

CPS Cordless Prosthodontic Screwdriver



Turns the implant world

save time | earn precision

bredent
medical

CPS Cordless Prosthodontic Screwdriver

Perfect Handling

- cordless
- ergonomic design

Highest precision


- Precision of torque values +/- 5 %
- Torque range: 8 Ncm – 40 Ncm

Best hygiene

- Lubefree handpiece and removable switch
- autoclavable
 - thermo-disinfectable

Universal

- With screw drivers for the most important implant systems



How much time do you spent loosening and fixing implant screws?

You can save app. 50 % of the time using the CPS.

Are you sure that your screws are always tight correctly?

With CPS there will be no screw loosening due to wrong torque.

Can you reach your palatinal screws easily?

With CPS you have a perfect cordless access to all sites, because it is cordless.

Order information:

CPS Cordless Prosthodontic Screwdriver

- Cordless motor
- Battery runtime for 40 screws
- Charging point
- Handpiece 80:1
- Screw drivers:

Hexagon

0,05 0,9 1.0 1.2 1.8 2.5

Torx

SKY 5.5 6.0

Slot

1.6 2.0

REF 580 CPS4 0

